Dietary fibre Filtration Module

(Ref. No: NIN/ST/12/Dietary Fibre/2017-18/7300-7311 &7398)

Filtration and dehydration capacity 6 samples/batch.

Filtration rate: 2 – 4 minutes/100 ml of digested solution,

depending on sample type.

Repeatability SD <0.5% absolute

Flask volume 600 ml

Accuracy According to official methods

Water Bath

Capacity 12 flasks

Temperature range Ambient to 80 °C

Temperature accuracy ±0.1°C

Shaking frequency 0 - 150 strokes/minute

Stroke length 0-50 mmPower supply 220 V,50Hz